



Certificate of Analysis

CP-Buffer-22

AxyGen's ready to use buffers, precipitants, salts, additives, and detergents provide the most complete and economical solution to crystallization of proteins, peptides, small molecules and nucleic acids. They can be used as a rapid Initial screening method or to accelerate the creation of optimization/follow up screens. Reagents are supplied in sterile bottles (200ml/bottle).

Product Name: 1.0 M MOPS, pH 6.5
Description: MOPS titrated with NaOH
Catalog #: CP-Buffer-22
Chemical Name: 3-(N-Morpholino)propanesulfonic acid
Molecular Formula: C₇H₁₅NO₄S
Molecular Weight: 209.26
CAS Number: [1132-61-2](#)
Beilstein Registry: 1106776
EC Number: 2144785
Purity: ≥99.5%
Appearance: White crystalline powder
Solubility: This product is soluble at 1g/2ml of water, yielding a clear, colorless solution
pKa (25 °C): 7.2



AxyGen BioSciences
33210 Central Ave., Union City, CA 94587
Ph.: (510) 476-1090 Fax: (510) 476-1091
Email: sales@axygenbio.com
Website: www.axygenbio.com